











Connector Guide



Manufacturers Models	EnviteC Part Number	Finger Clip Sensor OEM Part Number	Sensor Cable Length	Sensor Conn. No.	Adapter Cable Part Number	Adapter Cable Length	Extension Cable Part Number	Extension Cable Length	Adapt./Ext. Cable Conn. No. (Monitor side)	Adapt./Ext. Cable Conn. No. (Probe side)
Agilent Technologies (see Philips Medical)										
Armstrong										
AD-900, 1000, 2000, 3000, 4000	F-3222-12	BCI® 3044	1.20	C 03	direct connect					
Artema										
MM200	F-2712	[not available]	3.00	C 12	direct connect					
Diascope	F-3212-9	Nellcor® DS100A	0.90	C 05	5613	3.00			C 16	C 60
Diascope	F-2613		3.00	C 16	direct connect					
S&W	F-2612	S&W® SL033057	3.60	C 15	direct connect					
Baxter										
ASAT	F-3812	SiMED® 600015	3.60	C 23	direct connect					
Explorer/Monarch	F-3212-9	Nellcor® DS100A	0.90	C 05	direct connect		4212	2.40	C 05	C 60
BCI (Smiths Medical)										
3100, 3101, 6100, 9100, Advisor, Autocorr, Mini Torr, Mini Torr Plus	F-3222-12	BCI® 3044	1.20	C 03	4223	2.40			C 33	C 60
same as above	F-3222-30		3.00	C 33	direct connect					
3300, 3301, 3302, 3303, 3304, 3401	F-3222-12	BCI® 3044	1.20	C 03	direct connect		4222	2.40	C 03	C 60
Biosys										
BPM-200	F-3212-9	Nellcor® DS100A	0.90	C 05	direct connect		4212	2.40	C 05	C 60
Bird										
4400	F-3282-36	CSI® 511-10L	3.60	C 31	direct connect					
same as above	F-3282-9	CSI® 511SD	0.90	C 08	4282-1	2.40	4282	2.40	C 31	C 60
Bruker (see Odam)										
Burdick (see Sacelabs Burdick)										
CAS Medical Systems										
OscilloMate 9002, 9303, 9302S	F-3012	Nonin® 8000AA-1	0.90	C 04	direct connect		4012-1	3.00	C 04	C 60
OscilloMate 9002, 9303, 9302S	F-3012	Nonin® 8000AA-1	0.90	C 04	direct connect		4012	0.90	C 04	C 60
Colin										
Press-Mate SAT, BP 88, BP-306, BP-508, BP-808, BP-88S	F-3212-9	Nellcor® DS100A	0.90	C 05	4213	3.00	4212	2.40	C 05	C 60
Corometrics (see GE Medical Systems)										
CSI (Criticare System)										
501, 503, 504, 507, 508, 581, 1100, 8100, Poet (602-3A), Scholar	F-3282-9	CSI® 511SD	0.90	C 08	4282-1	2.40			C 31	C 60
(same as above)	F-3282-30		3.00	C 08	4282-1	2.40			C 31	C 60
(same as above)	F-3282-36	CSI® 511-10L	3.60	C 31	direct connect					
503 Spot, 504Plus, 504DX, 506DXN	F-3282-9	CSI® 511SD	0.90	C 08	direct connect		4282	2.40	C 08	C 60



503 Spot, 504Plus, 504DX, 506DXN Critikon (see GE Medical Systems) Datascope
Accusat, Accutorr 3/4 SAT, Accutorr Plus, (Datascope SpO ₂), Passport (Datascope SpO ₂) (same as above)
Expert, Passport (Nellcor® SpO ₂) Accutorr Plus (Nellcor® SpO ₂) Datex-Ohmeda
AS/3; CS/3; CH/S; CH-RS; CG/2S; Cardiocap I, II, 5; Capnomac Ultima; Light Monitor; M-Estpr; Oscar II; OSP-200; Satellite Trans Plus; ULTSV I; ULTSVO; ULTSIO; ULTS00 (same as above)
Ohmeda 3700, 3710, 3770, 3775; Biox 3740, 3760; 3800, 3900, 4700, TuffSat, OxiCap, 5250RGM, Rascal II (same as above)
Digicare Digimax 5500 Dima Italia All Dräger Medical Dräger Oxipac 2000 Vitalert 1000, Cicero Narkomed Elmed PSO1000 Entri Medical Propaq 102, 104, 106 EnviteC PulsiQuant 8040 Fukuda Denshi HG302 Module, DS-5100E, DS-5300W Gambro BP806

F-3282-30		3,00
F-3312-9	Datascope® 0600-00-0026-01	0.90
F3312-36	Datascope® 0020-00-071-01	3.60
F-3212-9	Nellcor® DS100A	0.90
F-3212-9	Nellcor® DS100A	0.90
F-3512-9	Datex® SAS-F	0.90
F-3512-40	Datex® SAS-F4	3.60
F-3412	Ohmeda® 6051-0000-035	3.60
F-3412-9	[none offered]	0.90
F-3222-12	BCI® 3044	1.20
F-3212-9	Nellcor® DS100A	0.90
F-3222-12	BCI® 3044	1.20
F-3212-9	Nellcor® DS100A	0.90
F-3212-9	Nellcor® DS100A	0.90
F-3222-12	BCI® 3044	1.20
F-3212-9	Nellcor® DS100A	0.90
F-3222-12	BCI® 3044	1.20
F-3212-9	Nellcor® DS100A	0.90
F-3212-9	Nellcor® DS100A	0.90

C 08	direct connect		4282	2.40	C 08	C 60
C 09	4312-1	3.00			C 51	C60
C 51	direct connect					
C 05	direct connect		4212	2.40	C 05	C 60
C 05	4213	3.00			C 41	C 60
C 01	4512	3.00			C 14	C 60
C 14	direct connect					
C 22	direct connect					
C 01	2241-8	2.40			C 22	C 60
C 03	4222-15	1.50	4222	2.40	C 33	C 60
C 05	direct connect		4212	2.40	C 05	C 60
C 03	direct connect		4222	2.40	C 03	C 60
C 05	5712	2.40			C 39	C 60
C 05	4214	3.10			C 38	C 60
C 03	direct connect		4222	2.40	C 03	C 60
C 05	4214	3.10			C 38	C 60
C 03	direct connect		4222	2.40	C 05	C 60
C 05	direct connect		4212	2.40	C 05	C 60
C 05	direct connect					

SpO₂ Cross Reference Guide




Connector Guide




Manufacturers Models

Manufacturers Models	EnviteC Part Number	Finger Clip Sensor OEM Part Number	Sensor Cable Length	Sensor Conn. No.	Adapter Cable Part Number	Adapter Cable Length	Extension Cable Part Number	Extension Cable Length	Adapt./Ext. Cable Conn. No. (Monitor side)	Adapt./Ext. Cable Conn. No. (Probe side)
GE Medical Systems										
7000 Series, TRAM A.R., TRAM 200SL, 400SL, 600SL, 800A, 800SL	F-3412	Ohmeda® 6051-0000-035	3.60	C 22	direct connect					
(same as above)	F-3412-9	[none offered]	0.90	C 01	2241-8	2.40			C 22	C 60
Eagle, Dash, SOLAR, TRAM Modular 250SL, 450SL, 650SL, 850A, 850SL	F-3212-9	Nellcor® DS100A	0.90	C 05	5812	2.40			C 50	C 60
Hellige Eagle 1000, 3000, SMU, SMK, SMV	F-3212-9	Nellcor® DS100A	0.90	C 05	5813	3.00			C55	C 60
Corometrics 188, 129,129F, 128, 128F	F-3212-9	Nellcor® DS100A	0.90	C 05	5812	2.40	4212	2.40	C 05	C 60
Corometrics 120, 510, 511, 118	F-3212-9	Nellcor® DS100A	0.90	C 05	5912	1.80			C 27	C 60
Critikon Dinamap Oxytrak, 8280	F-3112	Critikon® 8997	0.90	C 02	[none offered]					
Critikon Dinamap 8700 / 8710 / 8720 Plus	F-3912-1	Critikon® 9084 + 8731	3.00	C 32	direct connect					
Critikon Dinamap 8700 / 8710 / 8720 Plus, Oxyshuttle	F-3912	Critikon® 9084	0.90	C 46	4912	3.00			C 32	C 61
Critikon Dinamap 9700 / 9710 / 9720 Plus, Compact, MPS, Select	F-3212-9	Nellcor® DS100A	0.90	C 05	direct connect		4212	2.40	C 05	C 60
Dinamap Pro 100-400	F-3212-9	Nellcor® DS100A	0.90	C 05	4213	3.00	4212	2.40	C 05	C 60
Goldway Industrial										
UT4000A	F-3222-12	BCI® 3044	1.20	C 03	4222-15	1.50			C33	C 60
Hellige (see GE Medical Systems)										
Healthdyne										
All	F-3012	Nonin® 8000AA-1	0.90	C 04	direct connect		4012-1	3.00	C 04	C 60
Hewlett Packard (see Philips Medical)										
Innomed										
Oxycard, Innocare-T0	F-3012	Nonin® 8000AA-1	0.90	C 04	direct connect		4012-1	3.00	C 04	C 60
Invacare										
All	F-3222-12	BCI® 3044	1.20	C 03	direct connect		4222	2.40	C 03	C 60
Invivo Research										
4500, 4500 Plus, Millenia, Omega 1445, Omni-Trak TVS, Sat Scout	F-2912	Invivo® 9383	3.00	C 32	direct connect					
4500 Plus 2	F-3212-9	Nellcor® DS100A	0.90	C 05	4412	2.40			C 32	C 60
Ivy Biomedical										
405A, 405D	F-3212-9	Nellcor® DS100A	0.90	C 05	4214	3.10			C 38	C 60
405P, 405T	F-3212-9	Nellcor® DS100A	0.90	C 05	direct connect		4212	2.40	C 05	C 60




C31




C32




C33



C34



C35




C36




C37



C38



C39



C41

Kontron
7840, 7845, Minimom 7137, Micromon, Supermon, Kolormon, Trakmon
Lynn Medical
All
Magna Medical
MON-MTEK-30, 33, 35, 40, 85, 1200, 1201, 1300, 1401, 1406, 1411, 1413, 1501, 1503, MON-MTEK-MOD2, -MOD3
Mallinckrodt
N-100 (w / High Sat), N-200
(same as above)
NPB-290/5, N.390, N.3000, N-6000 (w / Oxismart)
(same as above)
All (not listed above)
All (not listed above)
Marquette (see GE Medical Systems)
Mediana
YM-2002 Series
MEDLAB
NANOX2
POX 10
Medsource
BPM-200
M.D.E.
Escort 100 OPT11A, Escort 300A
Escort II OPT30; Prism
MEK
MP-100, MP-110, MP-400, MP-1000
Mennen
Horizon XL (Ohmeda SpO ₂)
(same as above)
Horizon XL (Nellcor SpO ₂), Mercury 1200, MR 1303
Nellcor (see Mallinckrodt)







F-2812	Kontron® 0608010	3.60
F-3222-12	BCI® 3044	1.20
F-3222-12	BCI® 3044	1.20
F-3212-9	Nellcor® DS100A	0.90
F-3212-31	Nellcor® DS100A + EC8	3.10
F-3212-9	Nellcor® DS100A	0.90
F-3212-31	Nellcor® DS100A + EC8	3.10
F-3212-9	Nellcor® DS100A	0.90
F-3212-31	Nellcor® DS100A + EC8	3.10
F-3212-9	Nellcor® DS100A	0.90
F-3212-9	Nellcor® DS100A	0.90
F-3222-12	BCI® 3044	1.20
F-3222-36		3.60
F-3212-9	Nellcor® DS100A	0.90
F-3212-9	Nellcor® DS100A	0.90
F-2012		3.00
F-3412	Mennen® 0800-060-540	3.60
F-3412-9	[none offered]	0.90
F-3212-9	Nellcor® DS100A	0.90

C 11	direct connect					
C 03	direct connect		4222	2.40	C 03	C 60
C 03	direct connect		4222	2.40	C 03	C 60
C 05	4214	3.10			C 38	C 60
C 05	4214	3.10			C 38	C 60
C 05	4213	3.00			C 41	C 60
C 05	4213	3.00			C 41	C 60
C 05	direct connect		4212	2.40	C 05	C 60
C 05	direct connect		4212	2.40	C 05	C 60
C 05	DOC-10				C 41	C 60
C 03	direct connect		4222	2.40	C 03	C 60
C 34	direct connect					
C 05	direct connect		4212	2.40	C 05	C 60
C 05	4214	3.10				
C 05	direct connect		4212	2.40	C 05	C 60
C 06	direct connect					
C 22	direct connect					
C 01	2241-8	2.40			C 22	C 60
C 05	4214	3.10			C38	C 60

SpO₂ Cross Reference Guide



Connector Guide



Manufacturers Models

Nihon Kohden
AL-800PA Module; BSM-1101, 2101/2, 3101, 7100, 7200; OLV-1100A, Lifescope
(same as above)
AL-900P Module; BSM-1102, 9510; OGS-2000
BSM-4104A, 4114A
BSM-2301K, 2303K, 4103, 4113, 9510J/K
Nonin
8500
8600, 8700, 8800
Novamatrix (Respironics)
510, 511, 505, 515B/C, 520, 7000, 7100, Oxypleth
(same as above)
500
Handheld 512, 513, 710, 715,
Odam
All
Ohmeda (see Datex-Ohmeda)
Pace Tech
Vitalmax 800 Series (black face, Nellcor® SpO ₂)
Modell 520
Palco
All
All
Philips Medical
78352C, 78354C, 78834C, M1020A, M1025B, M1025A, M1350B, M1722A/B, M1732A/B, M2475B
(same as above), ViridaA1, A3
(same as above)
Physio-Control
Lifepak 12

EnviteC Part Number	Finger Clip Sensor OEM Part Number	Sensor Cable Length
F-2202-1	N.K.*TL-101S	3.00
F-2202-9	[none offered]	0.90
F-3212-9	Nellcor® DS100A	0.90
F-3212-9	Nellcor® DS100A	0.90
F-2203-16	N.K.*TL-101T	1,60
F-3012	Nonin® 8000AA-1	0.90
F-3012	Nonin® 8000AA-1	0.90
F-3262	Novamatrix® 8776-00	3.00
F-3262-9	[none offered]	0.90
F-3262-1	Novamatrix® 8660-00	3.00
F-3263		0.90
F-3222-12	BCI® 3044	1.20
F-3212-9	Nellcor® DS100A	0.90
F-3712	PaceTech® 4510 + 4530	3.60
F-2112	[not available]	0.90
F-2112-24	[not available]	2.40
F-2412	Nellcor® DS100A + Agilent® M1900B	3.60
F-3212-9	Nellcor® DS100A	0.90
F-2414-15	M1191A	1.50
F-3212-9	Nellcor® DS100A	0.90

Sensor Conn. No.	Adapter Cable Part Number	Adapter Cable Length	Extension Cable Part Number	Extension Cable Length	Adapt./Ext. Cable Conn. No. (Monitor side)	Adapt./Ext. Cable Conn. No. (Probe side)
C 18	direct connect					
C 01	2201-8	2.40			C 18	C 60
C 05	direct connect		4212	2.40	C 05	C 60
C 05	4213	3.00			C 41	C 60
C 10	[none offered]					
C 04	direct connect		4012-1	3.00	C 04	C 60
C 04	direct connect		4012	0.90	C04	C 60
C 21	direct connect		4613	3.00	C 21	C 26
C 01	2261-8	2.40			C 21	C 60
C 21	direct connect		4613	3.00	C 21	C 26
C 07	direct connect					
C 03	5112	3.00			C 30	C 60
C 05	4214	3.10			C 47	C 60
C 47	direct connect					
C 43	direct connect					
C 43	direct connect					
C 17	direct connect					
C 05	5412	3.00			C 17	C 60
C 19	5414	2.00			C 17	C 62
C 05	direct connect					

